




# High Performance Liquid Chromatography

## Sample Information

Sample Name:	BVT 2733	Injection Volume:	2µl
Operator:	Jefferson	Run Time:	20min
Injection Date:	5/23/2019 4:07:30 PM	Channel:	280nm
Batch No.:	GC157381	Acq Method:	T

## Result



Rank	RetTime(min)	Area(mAU*s)	Height	Area %
1	3.678	14.64661	1.50954	0.1633
2	7.634	8853.85059	1479.49.023	98.6932
3	7.971	29.818.53	7.01827	0.3324
4	8.229	20.63570	5.45152	0.2300
5	10.367	37.21034	9.04578	0.4148
6	12.519	14.920.97	2.90615	0.1663