




# High Performance Liquid Chromatography

## Sample Information

Sample Name:	Bz-Nle-Lys-Arg-Arg-AMC	Injection Volume:	10µl
Operator:	Jefferson	Run Time:	20min
Injection Date:	9/23/2019 4:07:12 PM	Channel:	248nm
Batch No.:	GA212211	Acq Method:	T

## Result



Rank	RetTime(min)	Area(mAU*s)	Height	Area %
1	8.201	43226	6559	0.52
2	9.704	22136	3827	0.27
3	10.569	32275	4831	0.39
4	12.130	11990	2166	0.15
5	12.349	8150602	1022006	98.67